

# Preventing Biofilm Formation on Implants: The *In Vitro* Performance of Surgical Irrigation

Fahimeh Tabatabaei, PhD, Rebecca McMahon, PhD, Nina Bionda, PhD

## INTRODUCTION

Surgical site infections (SSIs) and periprosthetic joint infections (PJIs) are serious complications that can arise following arthroplasty. Infection remains the leading cause of revision surgeries, and with the increasing number of total joint arthroplasty (TJA) procedures year over year, effective infection control has never been more critical.<sup>1,2</sup>

One solution to reducing the risk of these complications is the use of surgical irrigation, designed to remove microbial contamination from the tissues and thus prevent colonization of implantable materials and subsequent potential for biofilm formation. Biofilms are comprised of aggregates of microorganisms and extracellular polymeric matrix, they adhere to implant surfaces, are notoriously difficult to treat with antibiotics, and significantly increase the risk of PJIs. A growing body of evidence suggests that optimizing irrigation solutions could be a game changer in reducing infection rates in joint replacement surgeries.<sup>3,4</sup>

This white paper presents a comparative *in vitro* analysis of the biofilm prevention effectiveness of nine commercially available irrigation solutions on four commonly used implantable materials. These solutions utilize a range of antimicrobial technologies, including acids, surfactants, polyhexamethylene biguanide (PHMB), povidone-iodine (PVP-I), chlorhexidine gluconate (CHG),


benzalkonium chloride (BZK), and a triple antibiotic solution (TAS). Each solution was evaluated using a two-minute pre-treatment contact time on four different implantable materials. Two common bacteria, Gram-positive *Staphylococcus aureus* and Gram-negative *Pseudomonas aeruginosa*, known to cause surgical site infections, were used in simulated synovial fluid to mimic the conditions typically encountered in joint replacement surgeries.

## METHODOLOGY

Round 0.5 inch diameter coupons of four different implantable materials (BioSurface Technologies Corporation, Bozeman, MT) without surface modifications were used in this testing: stainless steel (SS), cobalt-chromium alloy (CoCr), titanium alloy (Ti-6AL-4V), and ultra-high molecular weight (UHMW) polyethylene. Two pathogenic bacteria were used as model organisms, Gram-positive methicillin-resistant *Staphylococcus aureus* (MRSA) ATCC USA300 and Gram-negative *Pseudomonas aeruginosa* ATCC 27312. The simulated synovial fluid inoculation medium was purchased from Biochemazone (item# BZ183 without preservatives). **Table 1** shows the irrigation solutions evaluated in the study.

To evaluate the effectiveness of commercial surgical irrigation solutions in preventing attachment of bacteria, coupons (n=4 per group), attached to the bottom of 24-well plates using

a medical grade silicone adhesive, were first exposed to surgical solutions for two minutes. After aspiration,  $1 \times 10^6$  CFU of bacteria in synovial fluid was added to each coupon. The inoculated coupons were incubated at 37°C for 24 hours. Following 24 hours of contact time with the inoculum, the coupons were detached from the plate, and bacteria adhered to the material were recovered through a series of sonication debridement steps. The bacteria were enumerated using standard serial dilution and agar plating method. A schematic of this method is shown in **Figure 1**.



**Figure 1:** Schematic of the biofilm prevention model used in this study.


| SOLUTION                                   | KEY INGREDIENTS                                                           | NOTES                                                                                                                                |
|--------------------------------------------|---------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Saline                                     | 0.9% sodium chloride Solution                                             | Commonly used in wound care to wash and cleanse tissues without antimicrobial activity                                               |
| Surgiphor™ Antimicrobial Irrigation System | 0.5% povidone-iodine (PVP-I)                                              | Potential cytotoxic effects of iodine on human cells; deep brown color and strong smell                                              |
| BIASURGE® Advanced Surgical Solution       | 0.1% polyhexamethylene biguanide (PHMB), chelators (EDTAs), vicinal diols | Synergistic composition with enhanced antimicrobial properties of a broad-spectrum, safe and effective antimicrobial compound        |
| Prontosan® Wound Irrigation Solution       | 0.1% polyhexamethylene biguanide (PHMB)                                   | Safe and effective antimicrobial compound                                                                                            |
| Irrisept® Antimicrobial Wound Lavage       | 0.05% chlorhexidine gluconate (CHG)                                       | Significant and irreversible damage to cell; cause tissue irritation; can cause allergic reactions                                   |
| Bactisure® Wound Lavage                    | 0.13% benzalkonium chloride (BZK), 10% ethanol, 0.6% acetic acid          | High cell toxicity even after brief exposure (1 minute) in a study on 3D cultures of human fibroblasts                               |
| Vashe® Wound Solution                      | 0.033% hypochlorous acid (HOCl)                                           | Significant erosion and wear on cobalt chrome and titanium surgical implants; bleach smell, limited efficacy in established biofilms |
| Xperience™ Advanced Surgical Irrigation    | 32.5 g/L citric acid, 31.3 g/L sodium citrate                             | May cause irritation or damage to sensitive tissues, reduced efficacy in established infections                                      |
| Triple Antibiotic Solution (TAS)           | gentamycin, bacitracin, cefazolin                                         | Risk of antibiotic resistance development, and it may be less effective against biofilm in short contact times                       |

**Table 1:** Products tested in this study and their ingredients


RESULTS

The results of the assays are shown in **Figures 2 and 3**. BIASURGE Advanced Surgical Solution (BIASURGE) demonstrated superior and consistent reductions in biofilm formation for both *Staphylococcus aureus* and *Pseudomonas aeruginosa* across all tested implantable materials. This remarkable performance, with over a 6-log reduction in microbial load, especially

against challenging pathogens like *P. aeruginosa*, underscores its role as a key component in infection prevention protocols designed to reduce costly revision surgeries. The simulated synovial fluid used in the study closely mimics clinical conditions encountered during joint replacement surgeries.<sup>5</sup> The fact that BIASURGE maintained its efficacy under these conditions further validates its use in real-world surgical applications.



**Figure 2:** Biofilm prevention effectiveness of various commercial surgical irrigation solutions on coupons after a twenty-four hour challenge with MRSA in simulated synovial fluid.



**Figure 3:** Biofilm prevention effectiveness of various commercial surgical irrigation solutions on coupons after a twenty-four hour challenge with *Pseudomonas aeruginosa* in simulated synovial fluid.

## DISCUSSION

Compared to the other PHMB-based solution (Prontosan), BIASURGE showed higher level of performance, which can be attributed to the synergistic antimicrobial formulation based on a proprietary blend of PHMB, EDTAs which serve as metal chelators, and vicinal diols which interact with lipids.<sup>6</sup>

The BZK/ethanol/acetic acid solution demonstrated favorable performance; however, there are ongoing concerns regarding its toxicity.<sup>7</sup> Additionally, while there is extensive corrosion testing data on metal-based materials for BIASURGE, comparable research on Bactisure is lacking and the recommended rinsing after application does not eliminate questions about its effects on metal implants during brief exposures.<sup>8</sup>

The material-dependent efficacy of Prontosan, Irrisept, and TAS underscores the complex interactions between irrigation solutions and specific implant materials. This complexity should be a key consideration when selecting a solution for clinical use.<sup>9</sup>

In contrast, the citric acid, HOCl, and PVP-I solutions were largely ineffective in preventing the formation of *S. aureus* and *P. aeruginosa* biofilms on any of the tested materials. While some studies have indicated that PVP-I solution is ineffective at eradicating *S. aureus*

biofilm, our findings expand on this by demonstrating that it also fails to prevent the formation of biofilms.<sup>10</sup> This dual ineffectiveness highlights the limitations of these solutions in clinical settings where biofilm-related infections are a significant concern, underscoring the need for more robust solutions like BIASURGE.

## CONCLUSION

The results of this study provide valuable insights into the effectiveness of various surgical irrigation solutions in preventing biofilm formation by *S. aureus* and *P. aeruginosa* on different implantable materials. One of the most significant findings is the superior performance of the PHMB-based technology in BIASURGE, which resulted in over 6-log reduction in microbial load across all implantable materials and pathogens tested compared to untreated control materials.

The study also highlights the importance of selecting an irrigation solution that is tailored to the specific implant materials used in surgeries. Solutions like BIASURGE, which consistently performed well across all materials, may be preferable for broad-spectrum protection. Conversely, other solutions, while effective in certain contexts, may require more strategic use depending on the specific implant material.

As antimicrobial resistance continues to rise and hospitals seek more cost-effective ways to reduce infections and improve patient outcomes, BIASURGE presents itself as a high-impact solution that can help mitigate these challenges while enhancing patient care. Further *in vivo* testing and clinical trials will help validate its long-term efficacy, thus improving patient outcomes in total joint arthroplasty and similar procedures.

## ACKNOWLEDGMENTS

This study was performed at iFyber in a randomized, blinded manner. We would also like to acknowledge Larry Estlack for his early contribution to data generation using a similar experimental model. His work laid an important foundation for the current study.

## RESEARCH REFERENCES

1. Eriksson HK, Lazarinis S. Patient-related factors associated with superficial surgical site infection and progression to a periprosthetic joint infection after elective primary total joint arthroplasty: a single centre, retrospective study in Sweden. *BMJ Open* 2022;12: e060754. doi:10.1136/bmjopen-2022-060754.
2. Brown ML, et al. Evolving etiologies, and rates of revision total knee arthroplasty: a 10-year institutional report. *Arthroplasty* (2022) 4:39. doi: 10.1186/s42836-022-00134-7.
3. Caid M, Valk J, Danoff J. Irrigation Solutions in Total Joint Arthroplasty. *Spartan Med Res J*. 2022;7(2):37502. doi: 10.51894/001c.37502.
4. Staats A, Li D, Sullivan AC, Stoodley P. Biofilm formation in periprosthetic joint infections. *Ann Jt*. 2021;6:43. doi: 10.21037/aoj-20-85.
5. Mian HM, Lyons JG, Perrin J, Froehle AW, Krishnamurthy AB. A review of current practices in periprosthetic joint infection debridement and revision arthroplasty. *Arthroplasty*. 2022;4(1):31. doi: 10.1186/s42836-022-00136-5.
6. Salamone AB, Salamone JC, McMahon RE, Poleon S, Bionda N, D'Arpa P. Synergistic Effect and Antibiofilm Activity of a Skin and Wound Cleanser. *Wounds*. 2020;32(8):208-216.
7. Sosnoski D, Dietz P, Bou-Akl T, Ren WP, Markel D. Irrigation Solutions Negatively Affect the Viability and Function of Human Fibroblasts: An in vitro Study. *Biomed Hub*. 2022;7(3):165-172. doi: 10.1159/000527110.
8. Rodier B, McMahon R, Bionda N. Compatibility Assessment of BIASURGE® Advanced Surgical Solution with Implantable Materials and Closure Technologies. *sanaramedtech.com*. May 2024.
9. Siddiqi A, Abdo ZE, Springer BD, Chen AF. Pursuit of the ideal antiseptic irrigation solution in the management of periprosthetic joint infections. *J Bone Jt Infect*. 2021;6(6):189-198. doi: 10.5194/jbji-6-189-2021.
10. Parker DM, Koch JA, Gish CG, Brothers KM, Li W, Gilbertie J, Rowe SE, Conlon BP, Byrapogu VKC, Urish KL. Hydrogen Peroxide, Povidone-Iodine and Chlorhexidine Fail to Eradicate *Staphylococcus aureus* Biofilm from Infected Implant Materials. *Life*. 2023; 13(6):1230. doi: 10.3390/life13061230